

AUXILIARY WALKING ELECTRIC WHEELCHAIR

MODEL NO.: **ZW518**




Free to stand from seat

Ergonomics design, press the rise-fall button to shift the sitting status to walking auxiliary status.



FEATURES

- Switch between wheelchair mode and standing mode with one button
- Help users stand and walk safely
- Assist patients with lower limb rehabilitation training
- Reduce the risk of bedsores caused by prolonged bed rest



One button to switch modes



Two rear drive motor



Electric-magnetic brake



Safety assurance in gait training

TECHNICAL SPECIFICATIONS

- Wheelchair mode size: 1000*680*1100mm
- Gait rehabilitation training mode size: 1000*680*2030mm
- Load capacity: 100KG
- Max. Drive Speed: 6 Km/h
- Maximum driving mileage with fully charged battery ≥ 20 km
- Tires Size: 8-inch front wheel, 10-inch rear wheel
- Applicable people: Height 140cm-190cm, Weight ≤ 100 KG
- Battery Charge Time: 4H
- Brake System: Electromagnetic brake
- Lithium Battery: Rated capacity: 15.6Ah; Rated voltage: 25.2V